Art. N° **63700**

ENG **ALCOHOL TEST**

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ENG ALCOHOL TEST

The alcohol test is intended for single use, it is used to test the alcohol content in the exhaled air. You can determine the alcohol content in the blood. The values on the comparative scale indicate: 0.0 ‰, 0.2 ‰, 0.5 ‰. The manufacturer is not legally responsible for the result or reliability of the test performed. The level of alcohol increases within an hour after its consumption. It should be remembered that even the smallest, permissible amount of alcohol in the blood may impair the ability to drive vehicles, thus posing a direct threat to the safety of life and health and road traffic hazards for its participants. The result of the alcohol test does not constitute an unambiguous confirmation of whether a person has the ability to drive a vehicle after drinking alcohol, without violating the law of offenses and criminal law in this regard. Alcohol test is intended for private use and is not a legalized product for measuring alcohol content. The alcohol test measurement cannot be used as evidence. Only the user of the alcohol test is responsible for the interpretation of the test result. The manufacturer is not responsible for using the product in a manner inconsistent with the operating instructions.

Instructions for use: In order to obtain the most accurate result, start measuring at the earliest 20 minutes after the last consumption of alcohol. Smoking, chewing gum, and the consumption of certain food and drink products may affect the accuracy of the measurement. Before using the alcohol test, it is not recommended to: wash and rinse the mouth with oral hygiene products.

1. Component parts. Alkoholtest consists of a glass tube, a measuring substance - a reagent, absorption granules, filters, aluminum foils and 2 protective caps. Important: Individual components of the alcohol test may be life or health-threatening. With improper use, the glass tube may break or break and cause cuts. After using the alcohol test, it should be thrown in the trash. Alcohol test does not pose a threat to the environment in accordance with the European Union Directive.

2. Usage. A.) The alcohol test should be opened immediately before use and should only be used once. B.) Press the protective caps on both sides of the alcohol test until the aluminum foil is pierced. C.) Take a deep breath and blow air through one side through the protective caps (air flow direction is not important) for at least 10 seconds. To be sure, you can blow in the air twice. D.) Wait at least two minutes until read the results. E.) If the color in the center of the tube changes from white to pink, the person is under the influence of alcohol. Compare the color intensity with the scale on the tube showing the amount of alcohol in the exhaled air. The most accurate comparison result can be obtained in daylight. F.) After 5 minutes of use - the result may be unreliable and the breathalyzer cannot be used again. G.) Dispose of the used alcohol test in the bin, it does not pose a threat to the environment in this form.

Attention! Taking into account the measurement error, the change in the color of the measuring substance after the test must not be ignored. The blowing of air cannot be interrupted or restarted! It should be remembered that any change in the comparative scale of the alcohol test means that the person taking the measurement is in a state of intoxication and must not drive a motor vehicle under any circumstances! In this case, in order to accurately measure the alcohol content in the blood, use a legalized breathalyzer.

Precautions

1. The alcohol test should be kept out of the reach of children.
2. It is forbidden to extract the measuring substance from the alcohol test.
3. It is recommended to use the alcohol test without first spraying the windshield of the car as the vapors may interfere with the result.
4. Avoid direct contact of the measuring substance with the eyes and skin.
5. Do not inhale or eat the measuring substance.
6. Do not throw the alcohol test into a fire.
7. If the measuring substance comes into contact with your eyes, flush your eyes with water as soon as possible and contact your doctor immediately.
8. In case of contact of the measuring substance with the skin, rinse the contact area with plenty of water. If you develop a skin reaction (redness, rash, allergic reactions, wounds) contact your doctor immediately.
9. In case of swallow the measuring substance, drink plenty of water, do not induce vomiting. Contact a doctor immediately.
10. In case of inhalation of the measuring substance, move to fresh air. Contact with doctor if necessary.

Nominal amount of the mixture 0.606 ml Potassium Iodate
CAS: 7758-05-6

Type of threat:

- H272 May intensify fire, oxidiser.
- H318 Causes serious eye damage.
- H315 Causes skin irritation.

Storage method

- Protect against sun and high temperature.
- Keep in a dry and safe place out of the reach of children.